Drägerwerk AG & Co.KGaA, 23542 Lübeck

Our reference
-
Phone number
+49 451 882-6166

Vulnerability and Incident Reporting
product-security@draeger.com

Dräger Product Security Advisory PSA-20-267-01

Windows Netlogon Vulnerability CVE-2020-1472 "Zerologon"

>> Dräger Medical Devices are not affected by "Zerologon".

To whom it may concern,

On September 18, 2020, the US department of homeland security warned in Emergency Directive 20-04[1] about security vulnerabilities present in several Microsoft Server operating systems.
These vulnerabilities, dubbed "zerologon", affect the Microsoft Windows Netlogon Remote Protocol (MSNRPC) and could allow an unauthenticated attacker with network access to a vulnerable Microsoft Windows Server with the domain controller role to completely compromise all Active Directory identity services.

None of Dräger medical devices or the accompanying infrastructure contains a Windows Server with the Domain Controller role.

[1] https://cyber.dhs.gov/ed/20-04/

The vulnerabilities can be fixed by applying the Microsoft-provided August security patches[2]. Affected operating systems can be patched normally, as soon as the monthly patch verification for the respective Dräger product has been published to ServiceConnect.

For reporting incidents and potential vulnerabilities in our devices, please refer to https://static.draeger.com/security to contact the Product Security team directly.

Best regards,

Dräger Product Security

Revisions

Version	Date	Change
1	2020-06-19	Initial Document

[2] https://portal.msrc.microsoft.com/en-US/security-guidance/advisory/CVE-2020-1472

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany
Postal address:
23542 Lübeck, Germany
Tel +49 451 882-0
Fax +49 451 882-2080
info@draeger.com
www.draeger.com
VAT no. DE135082211

Bank details:
Commerzbank AG, Lübeck
IBAN: DE95 2304 0022 0014 6795 00
Swift-Code: COBA DE FF 230
Sparkasse zu Lübeck
IBAN: DE15 2305 0101 0001 0711 17
Swift-Code: NOLADE21SPL

Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7903 HL
General partner: Drägerwerk Verwaltungs AG
Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7395 HL

Chairman of the Supervisory Board
for Drägerwerk AG & Co. KGaA
and Drägerwerk Verwaltungs AG:
Stefan Lauer
Executive Board:
Stefan Dräger (chairman)
Rainer Klug
Gert-Hartwig Lescow
Dr. Reiner Piske
Anton Schrofner